




DEEP BREATH

Contactless respiratory rate monitoring

Yekun Wang, Justin Trobec, Karen Wong
W251-001, Fall 2020

Respiratory rate is a sensitive signal for clinical worsening, but current assessments have downsides

- Manually counting breaths per minute
- Rough estimation
- Equipment: capnographers, oximeters



Contactless respiratory rate monitoring has several advantages

HOME




- Assist telemedicine assessments
- Use equipment patients already have at home

CLINICAL SETTING





- Reduce infectious or other hazardous exposure
- Reduce medical waste
- More natural breathing pattern

Monitoring using pose detection models can preserve privacy



No personal identifying data



$t_{001}, kp1...$

$t_{002}, kp1...$




$t_{003}, kp1...$

$t_{004}, kp1...$

...

$t_{450}, kp1...$

Goal: detect respiratory rate from keypoint motion




**16 breaths/min
(normal)**

Pose detection on Jetson with pre-trained ResNet18

- Simultaneous person detection and keypoint estimation
- Pre-trained on data from MSCOCO
- Our model: 8 upper body keypoints (eyes, ears, shoulders, nose, neck) and 3 distances (ear-shoulder x2, nose-neck)




Respiratory rate estimation using fast Fourier transform ensembled with CNN-based model



Training respiratory rate detection model



- >200 video clips annotated with respiratory rate
- Model trained on Jetson devices (NX, AGX) for 100 epochs

Model	Training MAE	Validation MAE
Frequency Estimation + Heuristics	n/a	10.89
Frequency Estimation + LR	8.44	7.91
FFT + Feedforward Network	6.31	9.50
CNN	3.42	7.37
Ensemble: FFT FNN + CNN	2.79	5.49



Results

- Validation MAE for 15-second clips ~5.49
- Example: 15-second clip of breathing at 15 BPM



Clip1	13.74
Clip2	15.48
Clip3	15.93
Clip4	19.38
Average	16.14

Demo



Breath Rate Detector
Live Streaming




Start Breath Rate Detection

Detector Status
Rate: 12.10 breaths/min

Edge Device (Xavier)

Docker Container

Web UI



Record

Respiration Rate Detector



Docker Container

mqtt

```
(topic: rate)
{
  subject:id,
  session:ts,
  interval:#,
  rate:#,
  frame1:img
}
```

Docker Container

forwarder

```
{
  subject:id,
  session:ts,
  interval:#,
  rate:#,
  frame1:img
}
```

AWS Instance (Amazon Linux)

Docker Container

mqtt

```
(topic: rate)
{
  subject:id,
  session:ts,
  interval:#,
  rate:#,
  frame1:img
}
```

Docker Container

Uploader

```
{
  subject:id,
  session:ts,
  interval:#,
  rate:#,
  frame1:path
}
```

AWS Instance (Amazon Linux)

Docker Container




Web UI

Patient: B324N
Date: 2020/11/14
Time: 08:00PM

12 B/min


DynamoDB

Tables:
sessions
subj_sessions



Implications

- Contactless monitoring can be cheap, accurate, and privacy-preserving
- Other applications: baby monitors, home sleep studies, seizure detection




Next steps

- Expand training examples
 - Different poses (side view, full body, lying down)
 - Include other types of motion (e.g., talking)
- Monitor multiple people in frame (e.g., disaster field hospitals)
- Optical flow for motion estimation may improve accuracy over keypoints alone

Questions?

EXTRA SLIDES



Respiration Rate Detector





$t_{001}, kp1, kp2, kp3 \dots kp9$
 $t_{002}, kp1, kp2, kp3 \dots kp9$
 \dots
 $t_{450}, kp1, kp2, kp3 \dots kp9$

Predict

14 breaths/min






14 breaths/min

Conv2D, 9 channels (keypoints)

-> Flatten

-> Dense

-> Output layer w/ linear activation



If keypoint not visible → fill with 0 or other constant?